



Using the UPWalker® Bridge

The UPWalker® Bridge is engineered to further enhance the stability of the larger UPWalker-L by coupling the two armrest supports together.

To couple or connect the Bridge when installed (see installation instructions below), simply pivot the left and right Bridge arms toward each other. Align the left arm so that it slips into the socket on the right arm. Ensure that the latch has clicked securely into place before using the UPWalker-L.


To uncouple or disconnect the Bridge, press the gray button in the middle of the Bridge latch. Hold the button down as you pivot the Bridge arms forward until they separate.

To install the Bridge, unscrew the gray knob from the U-shaped bracket on the end of each of the Bridge arms. Slide the U-shaped bracket over the front side of upright arm adjustment tube until it “pops” into place. Make sure the bracket is touching the bottom of the armrest when installed. Then reinsert the gray knob into the hole on the outside of the U-shaped bracket and turn clockwise to screw into the nut on the opposite side of the bracket until very tight and secure. Repeat for both Bridge arms. Check to confirm that the Bridge arms are fastened very securely on both tubes and don’t slide down the tubes. Verify that the left and right arms meet at the same height in the middle for proper coupling.



To remove the Bridge from your UPWalker-L:

1. Uncouple the two arms of the Bridge.
2. Unscrew the gray knob and remove the screw from the U-shaped bracket on the end of the Bridge arm.
3. Grip the Bridge arm close to the upright tube and pull off of the tube. You may gently rock the Bridge arm up and down to dislodge the arm from the tube.
4. Reinsert the gray knob into the hole on the outside of the U-shaped bracket and screw into the nut on the opposite side of the bracket until secure to avoid losing it.
5. Repeat these steps for both Bridge arms.



CAUTION: Do not put your full weight on the UPWalker Bridge or use it as a handle to carry the UPWalker. It is designed only for the purpose of connecting the two armrest supports, not for supporting the user or carrying the UPWalker.

Do not adjust the armrest height or attempt to fold the UPWalker-L with the Bridge arms closed.

When the Bridge is open the Bridge arms are more prone to damage. To avoid damage when the Bridge is open:

- Avoid pushing down or pulling up on the free end of the Bridge arms.
- Do not allow the UPWalker-L to fall over causing the free end of a Bridge arm to hit the ground or a hard object.

It is recommended that the Bridge be removed when the UPWalker-L is folded for storage or transport since the Bridge arms in the open position can be damaged more easily.

For customer support on the Bridge for your UPWalker-L, please contact LifeWalker Customer Service at:
(877) 488-0822 or info@lifewalkermobility.com