

## PainAway Pro with Heat

HCPCS:




FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Electrodes (available in SM and LG) provide mild electrical current for the temporary relief of pain from aching muscles
- Can be used for temporary pain relief in the shoulders, waist, back, neck, arms or legs
- 3 function options include TENS+heat, TENS only or Heat Only (2 temperature options) for customizable needs
- 8 modes with 20 intensity levels for user specific levels of stimulation
- Wireless device with remote control makes it easy to operate with hard to reach areas such as middle back.
- TENS unit contains rechargeable battery

PRODUCTS	Description
RTLAGF-1000	PainAway Pro with Heat 1/EA, 20/CV
SPECIFICATIONS	Description
Retail Box Dimensions	2 (L) X 6 (W) X 8 (H)
Retail Box Weight	0.75 lbs.
Warranty	Remote and Receiver--2 years from date of purchase. Gel electrode pads and cable are subject to normal wear and tear and are not covered under the warranty.
ACCESSORIES	Description
RTLAGF-1010	Electrode Pk 1LG/1SM, For RTLAGF-1000, PKG